

# Strategy for Diagnosis of Neuronal Ceroid Lipofuscinoses

Petra Oliva<sup>1</sup>, Thomas P. Mechtler<sup>1</sup>, Berthold Streubel<sup>2</sup>, David C. Kasper<sup>1</sup>

<sup>1</sup>ARCHIMEDlife GmbH, Vienna, Austria, <sup>2</sup>Medical University of Vienna, Department of Pathology, Vienna, Austria


## Introduction

The Neuronal Ceroid Lipofuscinoses (NCLs) are a group of inherited neurodegenerative disorders that affect children and adults. They are grouped together by similar clinical features and the accumulation of auto-fluorescent storage material. Enzyme deficiency NCLs are caused by genetic mutations resulting in enzyme deficiency leading to lipopigment accumulation and / or protein dysfunction necessary for cell mechanisms in nerve and other tissues. More than a dozen genes containing over 430 mutations have been identified to cause at least 13 known types of NCLs. The clinical differential diagnosis of the NCL types is based on age of onset, clinical phenotype, ultra-structural characterization of the storage material and enzyme levels. Symptoms associated with these disorders can vary widely. Although protein dysfunction or lipopigment accumulation influences many cells, brain cells are typically affected first. Clinical presentations include vision loss, epilepsy and myoclonic epilepsy, dementia, speech loss, movement disorder, behavior problems and learning disabilities and problems.

## Study Cohort

Our medical laboratory has developed, validated and accredited a novel diagnostic panel for differential diagnosis of NCLs utilizing a single Dried Blood Spot (DBS). Our assay includes testing for NCL1 and NCL2. Here we are presenting data from a high-risk population screening of over **1,185 cases suspicious of NCL1 and NCL2** by applying our developed diagnostic work-flow, including an MS/MS-based enzyme assay followed by genetic confirmatory testing using Next Generation Sequencing (NGS). We measured enzyme activity for palmitoyl protein thioesterase 1 (PPT1) and tripeptidyl-peptidase 1 (TPP1), followed by confirmatory genetic testing for over 60 cases.


Within this high-risk population screening, samples from **22 countries** were received, notably from Middle East (e.g. Saudi Arabia, Turkey) and Europe.




## Frequency of gene variants

| Gene | Nucleotide change   | Amino Acid change      | Frequency |
|------|---------------------|------------------------|-----------|
| PPT1 | c.271_287del17insTT | p.Gln91_Cys96delinsPhe | 1         |
|      | c.329A>G            | p.Asn110Ser            | 1         |
|      | c.451C>T            | p.Arg151*              | 1         |
|      | c.541G>T            | p.Val181Leu            | 1         |
|      | c.550G>A            | p.Glu184Lys            | 1         |
|      | c.559C>T            | p.His187Tyr            | 1         |
|      | c.677T>G            | p.Val226Gly            | 1         |
| TPP1 | c.1204G>T           | p.Glu402Ter            | 2         |
|      | c.1367C>A           | p.Ser456Tyr            | 1         |
|      | c.1417G>A           | p.Gly473Arg            | 1         |
|      | c.1547_1548delTT    | p.Phe516Ter            | 1         |
|      | c.225A>G            | c.225A>G               | 4         |
|      | c.379C>T            | p.Arg127Ter            | 3         |
|      | c.380G>A            | p.Arg127Gln            | 1         |
|      | c.466C>A            | p.Pro156Thr            | 1         |
|      | c.509-1G>C          | c.509-1G>C             | 4         |
|      | c.528delT           | p.Pro178GlnfsTer5      | 2         |
|      | c.622C>T            | p.Arg208Ter            | 3         |
|      | c.649G>C            | p.Gly217Arg            | 1         |
|      | c.697G>T            | p.Asp233Tyr            | 1         |

44 cases were submitted for **PPT1** genetic confirmation, and **8 patients were diagnosed**. 48 cases were submitted for **TPP1** genetic confirmation with final **diagnosis for 32 patients**.



The list of *PPT1* and *TPP1* variants identified is visualized below. The color and the size of the box represents frequency (bigger box represent more frequent variant).



## Conclusion

The presented data underlines the benefit of a fast and reliable diagnostic work-flow for NCLs in suspected individuals. The study was completed thanks to collaboration and support from BioMarin.